LUMENDI

DiLumen™

Endoluminal Interventional Platform
Plateforme interventionnelle endoluminale
Piattaforma per procedure endoluminali
Endoluminale interventionelle Plattform
Plataforma de intervención endoluminal
Plataforma intervencionista endoluminal

INSTRUCTIONS FOR USE
MODE D'EMPLOI
ISTRUZIONI PER L'USO
BEDIENUNGSANLEITUNG
INSTRUCCIONES DE USO
INSTRUÇÕES DE USO

LUMENDI

253 Post Road West

Westport, CT, 06880 • USA

info@lumendi.com • www.lumendi.com

1-844-481-0822

177-19751-00 REVISION H
07-2025

Emergo Europe B.V.
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

INTENDED USE

The DiLumen™ Endoluminal Interventional Platform (“DiLumen™” or “the device”) is an accessory to an endoscope. The DiLumen dual balloon accessory is intended for use with any standard endoscope. The device’s intended use is to ensure complete positioning of an endoscope during navigation in the large intestine, while assisting with optical visualization, diagnosis, tissue manipulation and endoscopic treatment.

SAFETY AND PRECAUTIONS

- Read and understand this manual before operating the device and save this manual for future reference. Failure to read and thoroughly understand the information presented in this manual may result in injury to the patient and/or user. In addition, failure to follow the instructions in this manual may result in damage to and/or malfunction of the device.
- DiLumen is a prescription medical device.
- Adult patients only (≥ 18 years old).
- Consult the medical literature relative to complications, hazards, and techniques prior to the performance of any endoscopic procedure.
- DiLumen™ should only be used by a physician trained in flexible endoscopy. A thorough understanding of the procedure, clinical applications and risks associated with flexible endoscopy is required before using the device in a procedure.
- It is the responsibility of each medical facility to ensure that only well-educated medical specialists, who have the proper training and competencies in endoscopic procedures and endoscopic equipment, will be involved in the use of DiLumen™.
- To ensure that the DiLumen can reach the target location, make sure the appropriate length device is used.
- Ensure that the outer diameter (“OD”) of the endoscope is compatible with the inner diameter (“ID”) of DiLumen™.
- Do not use the device if the package is open, damaged or expired. Use of a non-sealed or expired package might lead to injury, illness, or death of the patient.
- This device, including the inflation system, is a single-use device and was not designed to be reprocessed or re-used. LUMENDI cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety. Reprocessing or re-using the device may lead to device failure and/or contamination of the device and/or infection or cross-infection, which may result in injury, illness, or death of the patient.
- DiLumen™, including the inflation system, is not intended for re-use. Any institution, practitioner, or third party that reprocesses, refurbishes, remanufactures, sterilizes, and/or reuses this disposable medical device must bear full responsibility for its safety and effectiveness.
- Prior abdominal-pelvic surgery (i.e. hysterectomy, GU or colonic surgery), tortuous/floppy colon, hernias and stricture/tight angulation may make using the DiLumen more challenging.

- The balloons at the end of the disposable device can be punctured or damaged. Handle with care and avoid exposing them to sharp objects or rubbing them.
- All warnings, precautions and contraindications relevant to the endoscope are applicable to the use of the endoscope with the device. For any information required for the use of the endoscope, please refer to the corresponding instructions for use.
- Do not use this instrument for any purpose other than its intended use.

COMPLICATIONS

Possible complications that may result from the use of DiLumen™ include those relevant to standard endoscopy and other similar procedures. Possible complications include but may not be limited to perforation, hemorrhage, infection, or intestinal blockage.

CONTRAINDICATIONS

Contraindications include but are not limited to those specific to any endoscopic procedure, use of an overtube, or any endoscopic procedure performed in conjunction with an overtube, including: resistance in overtube advancement; colonic narrowing associated with but not limited to radiation injury, tumor, or inflammatory condition; angulated sigmoid colon; recent surgical alteration and/or anastomoses or rectal reconstruction; intractable patient discomfort; or inability to straighten the sigmoid colon. Not for pediatric use (< 18 years old).

CLINICAL BENEFITS

Clinical benefits include improved stability (and visualization) of the endoscope, reduction in colon looping via the balloons, as a conduit to efficiently remove, reinsert, or change endoscopes and the ability to manipulate tissue by extending/retracting the fore balloon during submucosal dissection via suture loops on the fore balloon.

ADVERSE AND SERIOUS INCIDENTS

The user, medical facility and / or patient should report all adverse events or serious incidents in relationship to the device to the Manufacturer, the device Importer / Distributor and the competent authority of the member State in which the user, medical facility and /or the patient is established.

Table of Contents

INTENDED USE	2
SAFETY AND PRECAUTIONS	2
COMPLICATIONS	3
CONTRAINDICATIONS	3
CLINICAL BENEFITS	3
ADVERSE AND SERIOUS INCIDENTS	3
I. INTRODUCTION	5
1. Purpose	5
2. Definitions and Glossary	5
3. Specifications	5
II. DEVICE DESCRIPTION	6
III. PREPARATION BEFORE THE PROCEDURE	8
IV. PERFORMING THE PROCEDURE	9
V. TISSUE MANIPULATION	13
VI. FLUSH PORT	14
VII. FINALIZING THE PROCEDURE	14
VIII. PLACING A SLIT IN THE SLEEVE	14
IX. DISPOSAL OF THE SINGLE-USE DILUMEN™ DEVICE	15
X. TROUBLESHOOTING	15
XI. STORAGE	16
XII. SYMBOL LEGEND	17

I. INTRODUCTION

1. Purpose

LUMENDI provides medical devices for endoluminal gastrointestinal procedures. This instruction manual describes how to operate the DiLumen™ device safely and effectively during a procedure within the large intestine. Before use, thoroughly review this manual. Keep this manual in a safe, accessible location.

2. Definitions and Glossary

Fore Balloon (FB): The front balloon of the device that can extend and retract in front of the endoscope.

Aft Balloon (AB): The rear balloon of the device that is fixed on the DiLumen™ sleeve.

Therapeutic Zone (TZ): The segment of the intestine that is isolated between the two inflated Fore and Aft balloons. The TZ length is the distance between the AB and the FB.

3. Specifications

- Storage environment: 0-40°C; Relative humidity: 15-85%
- Balloons (AB and FB) Outer diameter: 60 mm
 - Nominal pressure: 45 mmHg
 - Pressure range: 33-55 mmHg
 - Sleeve assembly: Inner diameter 17 mm


REF (qty per box)	Endoscope Working Length
D-1001 (-1,-5)	130 cm
D-1002 (-1,-5)	103 cm

- The DiLumen is a non-sterile device that does not require cleaning or sterilization prior to use.

II. DEVICE DESCRIPTION

The DiLumen™ device includes the following components:

1. Sleeve Assembly


- (1) Endoscope seal
- (2) Base
- (3) Fore Balloon (FB) knob
- (4) Sleeve

- (5) Push rods
- (6) Fore Balloon (FB)
- (7) Aft Balloon
- (8) Inflation lines

- (9) Inflation line connectors
- (10) Pinch Clamps
- (11) Suture Loops
- (12) Flush Port

The sleeve supports the two balloons (AB and FB) and is fixed onto the endoscope through the endoscope seal at the rear of the base. The FB is attached to two push rods, which can be manipulated (extended/retracted, locked/unlocked in place) using the FB knob on the base. The push rods are also used for inflation of the FB.

2. Inflation Handle


- (1) Selection switch
- (2) FB inflation indicator
- (3) AB inflation indicator
- (4) Inflation bulb
- (5) Inflation lines
- (6) Inflation line connectors


The inflation handle enables manual inflation or deflation of the two balloons independently by squeezing the inflation bulb. The inflation handle is connected to the sleeve assembly via the inflation lines and their connectors. The selection switch enables selection of each balloon and either inflation or deflation mode. The inflation indicators provide information on the inflation status of each of the balloons:

- Partially Red: deflated
- Mostly Yellow: partially inflated/deflated
- Only Green: fully inflated


Inflate AB


Deflate AB

Inflate FB

Deflate FB

NOTE

During the transition from deflated to inflated state or inflated to deflated state, the yellow color is visible and the indicator may oscillate before reaching a steady state

III. PREPARATION BEFORE THE PROCEDURE

The DiLumen™ device is a non-sterile, single use device that is supplied in a sealed tray. Before use, the device and packaging should be inspected to ensure no damage has occurred. If any damage is noted, **DO NOT USE**.

CAUTION

Do not use DiLumen™ if the package is open, damaged or expired. The use of a device from a compromised package may result in patient injury, illness, or damaged structural integrity of the device.

Be careful not to drop the sleeve assembly or the inflation handle (including inflation lines) when taking the device out of the package.

Connect the inflation lines securely to ensure proper inflation and deflation of the balloons.

Verify that nothing is excessively crushing, bending, or twisting the inflation lines to ensure proper inflation and deflation of the balloons.

1. Remove the sleeve assembly and inflation handle from the tray and uncoil the sleeve and the inflation lines. Lay the sleeve assembly straight and flat on a long plane surface.

2. Loosen the FB knob and fully advance the FB. Tighten the FB knob to lock the FB in place for insertion of the endoscope.

3. Lubricate the Sleeve. The inside of the DiLumen incorporates a self-lubricating coating. The coating is activated by hydrating with water or saline. The coating reduces friction between the endoscope and the sleeve lumen. To activate the coating execute the following:

Required materials:

- 30cc Luer-Lock™ dispensing syringe
- Sterile water or saline


Lubricating instructions:

- Fill the syringe with sterile water (min 30 cc).
- Hold the distal and proximal sleeve ends so the sleeve forms the shape of a U.
- Inject a minimum of 30cc of water or saline into the sleeve through the device handle.
- Move the distal and proximal sleeve ends up and down so the liquid contacts the sleeve's entire internal surface. Continue this hydrating process for a minimum of 30 seconds.
- After 30 seconds, discard the excess liquid through the distal end of the device.

4. Advance the endoscope through the sleeve until the AB is directly behind the articulation joint of the endoscope.

5. Make sure to orient the base so that the push rod control ramp is aligned with the endoscope tip movement wheels (in the same plane) as shown below.


6. Connect each of the two inflation lines from the inflation handle to the sleeve assembly using the connectors at the end of each line. The end of each line should be matched according to the shape of the connectors. The connection is unambiguous based on the male/female configuration.
7. Fully inflate and deflate each balloon once to make sure that the inflation system is working properly. Loosen the FB knob and partially extend and retract the FB to ensure that the FB may be manipulated properly. Do not fully retract the FB against the endoscope tip.
8. Apply water-based lubricant around the circumference of the tip of the endoscope and to the FB inner lumen. Fully retract the FB onto the tip of the endoscope. The endoscope tip should be approximately 10 cm in front of the FB. Tighten the FB knob to lock the FB in place on the endoscope shaft.
9. DiLumen™ is now ready for insertion.

IV. PERFORMING THE PROCEDURE

CAUTION	DiLumen™ should only be used by a physician trained in flexible endoscopy. Thorough understanding of the procedure, clinical applications and risks associated with flexible endoscopy is required before using the device in a procedure.
NOTE	If targeting multiple locations within the intestine, the most proximal location should be targeted first, followed by the second most proximal location, etc.

A – Navigation to target location

1. Verify that the two balloons are fully deflated.
2. Verify that the FB is docked on the endoscope tip and that the FB knob is tightened and fully locked.
3. Apply lubricant to the tip of the endoscope and to the outside of the FB and AB.
4. Insert the endoscope with the mounted device (FB and AB deflated) into the large intestine and navigate as performed in standard colonoscopy until the target location is reached.

CAUTION

Do not attempt to advance the endoscope and device without the FB fully docked onto the endoscope tip.

NOTE

Liberal use of insufflation and water irrigation is helpful when navigating the DiLumen/ Endoscope.

5. Connect each of the two inflation lines if not already connected.


B – Creation of the Therapeutic Zone (TZ)

CAUTION

The FB should always be advanced under direct visualization. If the TZ is deployed in a curved segment, see section B ii below for safe TZ deployment.

i. If the proximal side of the lesion is visible

6. Inflate AB to stabilize the endoscope by turning the selection switch to the *inflate AB* position and squeezing the bulb slowly until the AB inflation indicator is stable and shows green only. The AB Pinch Clamp may be closed at this time if desired.


NOTE

Squeezing the bulb slowly is more effective for inflating the balloons.

7. Verify proper AB anchoring by gently pulling on the scope and device and visually assessing grip of the colon through the endoscopic view.

8. Loosen the FB knob and extend the FB ahead of the endoscope using the FB knob until the desired distance between the FB and the endoscope tip is reached or until full extension of the FB is achieved.


CAUTION

The FB should be deflated when being advanced. Inflation of the FB before extension might result in damage to the device and/or inability to appropriately manipulate the FB.

If excessive resistance to the advancement of the FB is felt by the physician, the FB should not be pushed further into the intestine. The FB should be retracted and the endoscope articulation section should be straightened and reoriented before extending the FB again.

The FB should be retracted through the center of the balloon as damage to the pushrods may otherwise occur.

9. Inflate FB by turning the selection switch to the *inflate FB* position and squeezing the bulb slowly until the FB inflation indicator is stable and shows green only. The FB Pinch Clamp may be closed at this time if desired.

10. Verify proper FB anchoring by gently pulling on the FB knob and visually assessing grip of the colon wall through the endoscopic view. When satisfied with the position of the FB, tighten the FB knob to lock the FB in position.


11. The TZ is now established and the endoluminal procedure can be performed inside the TZ. Insufflate inside the TZ until a usual, acceptable endoscopic visualization is reached.

CAUTION


As in standard endoscopy, a minimum amount of gas should be used to maintain an acceptable field of view inside the intestine. An excess of insufflation in the TZ might result in over-inflation and injury to the patient.

ii. If the proximal side of the lesion is not visible (curved segment)


12. With both balloons deflated and the FB still docked on the endoscope tip, navigate the endoscope and device beyond the curved segment, proximal to (upstream of) the target lesion under direct visualization to a relatively straighter segment of the intestine.

13. Loosen the FB knob and extend the FB until it is entirely visible in front of the endoscope tip.


14. Inflate FB by turning the selection switch to the *inflate FB* position and squeezing the bulb slowly until the FB inflation indicator is stable and shows green only.


15. Verify proper FB anchoring by gently pulling on the FB knob and visually assessing grip of the colon wall through the endoscopic view.

16. With the FB now anchored, the FB knob loosened, and the AB still deflated, carefully withdraw the endoscope and device away from the FB, through the curved segment, and distal to (downstream of) the target lesion.

17. Tighten the FB knob to lock the FB in place. Inflate the AB to anchor and stabilize the endoscope by turning the selection switch to the *inflate AB* position and squeezing the bulb slowly until the AB inflation indicator is stable and shows green only.

18. The TZ is now established and the endoluminal procedure can be performed inside the TZ. Insufflate air into the TZ until a usual, acceptable endoscopic visualization is reached.

CAUTION

As in standard endoscopy, a minimum amount of gas should be used to maintain an acceptable field of view inside the intestine. An excess of insufflation used in the TZ might result in over-inflation and injury to the patient.

When using endoscopic tools within the TZ, take care to avoid contact with the FB. Sharp or hot objects can damage or puncture the FB and cause a loss of the TZ.

Do not attempt to retroflex an endoscope with the DiLumen™ device mounted. Use of the device affords the user the means to achieve superior endoscopic visualization.

C. Re-docking of the FB and deflation of AB

19. Once the endoluminal procedure is complete, open the FB Pinch Clamp if it was closed, deflate the FB by turning the selection switch to the *deflate FB* position and squeezing the bulb slowly under endoscopic visualization until the FB inflation indicator shows mostly red. The inner lumen of the FB should be clearly open and the FB visibly deflated.

20. Verify that the FB knob is adequately loosened. Align the tip of the endoscope with the FB inner lumen and re-dock the FB onto the tip of the endoscope by gently pulling the FB knob all the way back. With the selection switch turned to the *deflate FB* position, squeeze the bulb two times and verify that the FB inflation indicator is showing mostly red to ensure full deflation of the FB. Tighten the FB knob to lock the FB in place on the endoscope shaft.

CAUTION

The FB should be retracted through the center of the balloon as damage to the pushrods may otherwise occur.

If the FB cannot be successfully re-docked and/or the physician feels too much resistance during attempted FB re-docking, do not pull further on the FB knob. This may damage the push rod system.

Do not attempt to move the endoscope and the DiLumen device without fully opening both pinch clamps and deflating both balloons.

NOTE

If the FB cannot be successfully re-docked, the endoscope and device may be carefully withdrawn to a straighter segment of the intestine, which can facilitate easier re-docking. Deflate the AB and carefully withdraw the endoscope and device with direct endoscopic visualization of the un-docked FB.

If the FB still cannot be successfully re-docked, the FB may be overinflated or over-deflated. Ensure that the FB inflation level is adequate to allow for opening of the FB lumen and passage of the endoscope tip.

If the FB still cannot be successfully re-docked, the push rods may be crossed due to twisting of the device sleeve. Ensure that the sleeve is as straight as possible.

If the FB still cannot be successfully re-docked, the endoscope and device may be carefully removed with both balloons deflated under direct endoscopic vision.

21. Deflate the AB by opening the AB Pinch Clamp if it was closed, turning the selection switch to the *deflate AB* position and squeezing the bulb slowly until the indicator is showing mostly red for full deflation.

22. If needed, navigate the endoscope and device to the next most proximal target location and repeat these steps.

V. TISSUE MANIPULATION

1. The DiLumen balloons can be used to manipulate tissue inside the colon. With the AB inflated and anchored against the colon wall, the Fore Balloon can be extended and retracted as needed to reposition the colon, straighten Haustra folds, manipulate polyps, etc. to improve visualization and facilitate endoluminal interventions.

2. The Fore Balloon neck has two suture loops available to assist with tissue manipulation. For example, a loop can be used in combination with an endoscopic clip to attach a polyp to the loop. Once the polyp is attached to the loop, the push rods can be extended and retracted to pull the polyp for improved visualization or to retract the tissue during ESD or EMR.

CAUTION

The FB should always be extended and retracted slowly under direct visualization. When using the suture loops to manipulate tissue, the FB should be partially or fully deflated. Do not use excessive force to extend and retract the FB.

VI. FLUSH PORT

This DiLumen device includes a Flush Port in the Base. The Flush Port can be used to dispense sterile water or saline into the DiLumen Sleeve by connecting a filled luer-locking syringe to the Flush Port and pushing the syringe plunger. It is recommended that sterile water or sterile saline be used for flushing the DiLumen Sleeve.

VII. FINALIZING THE PROCEDURE

1. Verify that both balloons are deflated and carefully remove the endoscope and device from the patient.

NOTE

If during withdrawal of the endoscope the intestine wall needs to be examined or if an intervention is needed (e.g., biopsy, polypectomy, etc.), the two balloons or just the AB may be deployed again to stabilize the endoscope and facilitate the procedure. Deflate both balloons before further withdrawal of the endoscope.

2. Straighten the endoscope's tip and pull the endoscope out of the sleeve assembly completely.

CAUTION

Additional tools are not necessary in order to assist in loosening and removal of the device. Use of additional tools to facilitate removal may result in damage to the endoscope.

During the procedure, fluids might build up between the endoscope and the sleeve. Carefully remove the endoscope from the sleeve assembly to avoid fluids spillage.

VIII. PLACING A SLIT IN THE SLEEVE

The DiLumen Sleeve is a thin-walled thermoplastic polyurethane tube. If desired, once the DiLumen is at the target site, the inner lumen of the Sleeve can be accessed by using a scalpel and / or scissors to cut a slit in the Sleeve. This allows the user to create another access port into the Sleeve to use endoscopes or endoscopic tools during a procedure.

The slit should be made parallel to the length of the Sleeve and be centered between both Push Rods.

To form the slit, put a hole in the center of the Sleeve using a sharp instrument such as a scalpel or scissors close to the anus and then proceed to create a slit of about 2cm for access.

To prevent damage to the endoscope, it is recommended to remove the endoscope from the Sleeve prior to cutting the slit.

Use the smallest length slit necessary to minimize fluid and air leakage from the sleeve.

CAUTION

Be careful using sharp instruments such as scalpels or scissors when creating the slit to avoid injury to the patient or user.

Do not cut a slit in the sleeve until the endoscope and the Sleeve are at the target site to ensure the split is placed in the desired location on the sleeve. Verify that the slit does not have rough or sharp edges that could possibly cause tissue damage.

IX. DISPOSAL OF THE SINGLE-USE DILUMEN™ DEVICE

After use, the DiLumen™ single-use device may potentially be contaminated. Handle and dispose of the device (sleeve assembly and inflation handle) in accordance with accepted medical practice and applicable local, state and federal laws and regulations.


CAUTION

Do not re-use the single-use DiLumen™ device.

X. TROUBLESHOOTING

The Fore Balloon cannot be deflated

1. Check that the FB Pinch Clamp on the inflation line is not closed.
2. Check if the inflation lines between the base and the inflation handle are bent or kinked and if the inflation line connectors are properly connected. Attempt to deflate the FB by turning the selection switch to the *deflate FB* position and squeezing the bulb slowly until the FB inflation indicator is showing mostly red for full deflation. Assess FB deflation through the endoscopic view.
3. If the FB still does not deflate, de-couple the connectors and allow the FB to deflate. Assess FB deflation through the endoscopic view.
4. If the FB still does not deflate, make a hole in the FB under endoscopic view using a disposable endoscopic injection needle. Be careful not to cause tissue trauma. Assess FB deflation through the endoscopic view.
5. With both balloons deflated and any endoscopic tools removed from the endoscope, carefully remove the endoscope and device from the patient's body.

The Aft Balloon cannot be deflated

1. Check that the AB Pinch Clamp on the inflation line is not closed.
2. Check if the inflation lines between the base and the inflation handle are bent or kinked, and if the inflation line connectors are properly connected. Attempt to deflate the AB by turning the selection switch to the *deflate AB* position and squeezing the bulb slowly until the AB inflation indicator is showing mostly red for full deflation.
3. If the AB still does not deflate, deflate the FB by turning the selection switch to the *deflate FB* position and squeezing the bulb slowly until the FB inflation indicator is showing mostly red for full deflation. With the FB deflated, de-couple the connectors and allow the AB to deflate.
4. If the AB still does not deflate, carefully pass a small diameter endoscope alongside the existing endoscope and device. Once visualization of the AB is achieved, use an endoscopic tool to puncture the AB and allow deflation. Assess AB deflation through the endoscopic view of the small diameter endoscope.
5. Remove the small diameter endoscope. With both balloons deflated and any endoscopic tools removed from the endoscope, carefully remove the endoscope and device from the patient's body.

XI. STORAGE

The DiLumen™ device should be stored at room temperature under dry conditions.

CAUTION

Do not use the device after use-by date printed on label.

XII. SYMBOL LEGEND

	Medical Device
	Authorized Representative in the European Community / European Union
QTY	Quantity Symbol
	Federal (US) law restricts this device to sale by or on the order of a physician
	Consult Instructions for use or consult electronic instructions for use
	Use-by date
	Catalog Number
	Lot Number / Batch Code
	Do not use if package is damaged or open and consult instructions for use

	Do not re-use
	Temperature Limit
	Humidity Limitation
	Keep dry
	Medical Device Manufacturer
	Caution
	Non-sterile

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1. Verifique se o grampo de aperto AB na linha de inflação não está fechado.