LUMENDi

Pioneering Endosurgery.

DiLumen EZ¹

Endolumenal Interventional Platform
Plateforme interventionnelle endoluminale
Piattaforma per procedure endoluminali
Endoluminale interventionelle Plattform
Plataforma de intervención endoluminal
Plataforma intervencionista endoluminal

INSTRUCTIONS FOR USE (IFU) MODE D'EMPLOI ISTRUZIONI PER L'USO BEDIENUNGSANLEITUNG INSTRUCCIONES DE USO INSTRUÇÕES DE USO

 **LUMENDi**

253 Post Road West
Westport, CT, 06880 • USA
info@lumendi.com •
www.lumendi.com
1-844-481-0822

Emergo Europe B.V.
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

2797

265-06468 REVISION E
2025-04

THIS PAGE INTENTIONALLY LEFT BLANK

INTENDED USE

The DiLumen EZ™ Endoluminal Interventional Platform (“DiLumen EZ” or “the device”) is an endoscope accessory intended to ensure complete positioning of an endoscope in the large intestine, and assist with optical visualization, diagnosis, and endoscopic treatment.

SAFETY AND PRECAUTIONS

- Read and understand this manual before operating the device and save this manual for future reference.
- Failure to read and thoroughly understand the information presented in this manual may result in injury to the patient and/or user. In addition, failure to follow the instructions in this manual may result in damage to and/or malfunction of the device.
- Consult the medical literature relative to complications, hazards, and techniques prior to the performance of any endoscopic procedure.
- DiLumen EZ¹ is a prescription medical device.
- To be used in adult patients (≥ 18 years old).
- DiLumen EZ¹ should only be used by a physician trained in flexible endoscopy. A thorough understanding of the procedure, clinical applications and risks associated with flexible endoscopy is required before using the device in a procedure.
- It is the responsibility of every operator to ensure they have been properly trained prior to conducting any procedure using DiLumen EZ¹.
- To ensure that the DiLumen EZ¹ can reach the target location, make sure the appropriate length device is used. Options include 130 cm and 103 cm.
- Ensure that the outer diameter (“OD”) of the endoscope is compatible with the inner diameter (“ID”) of DiLumen EZ¹.
- Do not use the device if the package is open, damaged, or expired. Use of a non-sealed or expired package might lead to adverse health consequences.
- DiLumen EZ¹, including the Inflation System, is a single-use device and was not designed to be reprocessed or re-used. LUMENDI cannot verify that reprocessing can clean or maintain the structural integrity of the device to ensure patient and/or user safety. Reprocessing or re-using the device may lead to device failure and/or contamination and/or cross-infection, which may result in adverse health consequences.
- It is recommended to examine the patient’s large intestine with an endoscope alone prior to installing the DiLumen EZ¹ to determine if the anatomy presents any listed Warning, Precaution, Contraindication, or potential surgical challenges. Prior abdominal-pelvic surgery (i.e. hysterectomy, GU or colonic surgery), tortuous/floppy colon and hernias may make using the DiLumen EZ¹ more challenging.
- The balloon at the end of the disposable device can be punctured or damaged. Handle with care and avoid exposing the balloon to sharp objects or excessive manipulation.
- All warnings, precautions, and contraindications relevant to the endoscope are applicable to the use of the endoscope with the device. For any information required for the use of the endoscope, please refer to the corresponding instructions for use.
- Do not use this device for any purpose other than its intended use.

COMPLICATIONS

Possible complications that may result from the use of DiLumen EZ¹ include those relevant to standard endoscopy and other similar procedures. Possible complications include, but may not be limited to perforation, hemorrhage, infection, or intestinal blockage.

CONTRAINDICATIONS

Contraindications include, but are not limited to, those specific to any endoscopic procedure, use of an overtube, or any endoscopic procedure performed in conjunction with an overtube, including:

- Resistance in overtube advancement
- Colonic narrowing associated with but not limited to radiation injury, tumor, or inflammatory condition
- Angulated sigmoid colon
- Recent surgical alteration and/or anastomoses or rectal reconstruction
- Intractable patient discomfort
- Inability to straighten the sigmoid colon
- Not for pediatric use (<18 years old)

CLINICAL BENEFITS

Clinical benefits include improved stability (and visualization) of the endoscope, reduction in colon via the balloon and as a conduit to efficiently remove, reinsert, or change endoscopes.

ADVERSE AND SERIOUS INCIDENTS

The user, medical facility, and/or patient should report all adverse events or serious incidents related to the device to the Manufacturer, the device Importer/Distributor as well as the competent authority of the Member State in which the user, medical facility, and/or the patient is established.

Table of Contents

- 1. INTRODUCTION 6
- 2. DEVICE DESCRIPTION..... 7
- 3. PREPARATION BEFORE THE PROCEDURE..... 9
- 4. PERFORMING THE PROCEDURE..... 11
- 5. FINALIZING THE PROCEDURE..... 13
- 6. TROUBLESHOOTING 14
- 7. STORAGE..... 15
- 8. SYMBOL LEGEND 15

1. INTRODUCTION

1. Purpose

LUMENDI provides medical devices for endoluminal gastrointestinal procedures. This Instruction for Use (IFU) manual describes how to operate the DiLumen EZ¹ device safely and effectively during a procedure within the large intestine. Before use, thoroughly review this manual. Keep this manual in a safe, accessible location.

2. Definitions and Glossary

Inflation System: Comprised of the inflation handle, one balloon, inflation line and its connector, which are part of the Sleeve Assembly.

3. Specifications

- Storage environment:
 - 0-40⁰ C
 - Relative humidity: 15-85%
- Balloons:
 - Outer diameter: 60 mm inflated
 - Nominal pressure: 45 mmHg
 - Pressure range: 33 – 55 mmHg


- Product REF number:

REF (qty per box)	Endoscope Working Length
D-1100 (-1,-5)	130 cm
D-1101 (-1,-5)	103 cm

DiLumen EZ¹ is a non-sterile device that does not require cleaning or sterilization prior to use.

2. DEVICE DESCRIPTION


The DiLumen EZ¹ system consists of a **Sleeve Assembly (1)** and an **Inflation Handle System (2)**.

1. Sleeve Assembly Detail:


The Sleeve Assembly is secured onto the endoscope.

- (1) **Device Seal** – secures the sleeve to the endoscope
- (2) **Base** – supports attachment and stability
- (3) **Sleeve** – houses the endoscope
- (4) **Balloon** – inflates to stabilize positioning
- (5) **Inflation Line** – connects the balloon to the Inflation Handle
- (6) **Pinch Clamp** – locks inflation when engaged

2. Inflation Handle System Detail.

- The Inflation Handle System allows precise control of balloon inflation and deflation.
- The Balloon stabilizes the endoscope, assisting in positioning and procedural accuracy.

(1) **Connector** – links to the Sleeve


(2) **Inflation Line Connector**

(3) **Selection Switch** – toggles between inflation and deflation


(4) **Balloon Inflation Indicator** – shows inflation status
Red = Deflated
Yellow = Transitioning
Green = Fully Inflated

(5) **Inflation Bulb** – manually inflates and deflates the balloon

Inflate

Deflate

NOTE

When inflating or deflating, the indicator may briefly show yellow and fluctuate before stabilizing.

3. PREPARATION BEFORE THE PROCEDURE

CAUTION

Check the packaging: Do not use if the package is **open, damaged, or expired**. Using a compromised device may cause patient harm or device failure.

Handle with care: Keep the device clean and dry when removing it from packaging.

Verify Inflation Line Connection: Ensure the Inflation Line is securely connected to avoid inflation/deflation issues. Check that the Inflation Line is not bent, twisted, or crushed.

1. Uncoil & Lay Flat:

- Remove the **Sleeve Assembly** and **Inflation Handle** from the tray.
- Lay the **Sleeve and Inflation Line** flat.

2. Lubricate the Sleeve:

The DiLumen EZ1 has a self-lubricating coating that activates with water or saline.

Required Materials:


- 30cc Luer-Lock™ syringe
- Sterile water or saline

Lubrication Instructions:

- Fill the syringe with **at least 30cc of sterile water or saline**.
- Hold both ends of the **Sleeve in a U-shape**.
- Inject **30cc of liquid** through the device handle.
- Move both ends up and down for **30 seconds** to coat the interior.
- Drain any **excess liquid** through the distal end.

3. Insert the Endoscope:

- Insert a **clean endoscope** through the **Device Seal** into the Sleeve.
- Advance the endoscope until its tip aligns with the distal end of the Sleeve.

4. Connect the Inflation System:

- Attach the **Inflation Line Connector** to the **Sleeve Assembly**.

5. Test Balloon Inflation & Deflation:

- Fully **inflate and deflate the balloon** once to confirm functionality.

6. Final Lubrication:

- Apply a **water-based lubricant** to the **distal end of the Sleeve**.
- **Avoid excess lubrication** on the **Balloon**.

The device is now ready for use.

4. PERFORMING THE PROCEDURE

CAUTION

- Only physicians trained in flexible endoscopy should use this device.
- Ensure familiarity with the procedure, clinical applications, and associated risks before use.

NOTE

If targeting multiple locations within the intestine, the most proximal location should be targeted first, followed by the second most proximal location, etc.

A. Navigation to Target Location

1. Ensure the Balloon is fully deflated before insertion.
2. Lubricate the DiLumen EZ¹ tip before insertion.
3. With the Balloon deflated, insert the endoscope and navigate to the target location as in a standard colonoscopy.

NOTE


Use **ample insufflation and water irrigation** for smoother navigation.

B. Balloon Inflation

4. Position the target lesion in the desired monitor view (e.g., at 6:00).

5. To stabilize the endoscope, inflate the Balloon:

- Turn the **Selection Switch to inflate (+)**.
- **Squeeze the Inflation Bulb slowly** until the Inflation Indicator turns solid green.
- (Optional) Close the **Pinch Clamp** to maintain inflation.

NOTE

Slow, steady squeezing of the Inflation Bulb improves inflation control.

6. Check stability:

- Gently **tug the endoscope** to confirm anchoring.
 - **Visually assess** the grip through the endoscope
7. **Once stabilized, begin the procedure.** Use **minimal colon insufflation** to achieve acceptably clear visualization.

CAUTION

As is always the case with standard endoscopy procedures, avoid excessive insufflation to prevent over-inflation and patient injury.

Do not attempt to retroflex an endoscope while docked inside the Sleeve.


C. Ballon Deflation

CAUTION

Do **not** move the endoscope or DiLumen EZ¹ until the **Ballon is fully deflated and the Pinch Clamp is open.**

8. To deflate the Balloon:

- Open the **Pinch Clamp** (if closed).
- Turn the **Selection Switch** to **deflate (-)**.
- **Squeeze the Inflation Bulb slowly** until the **Inflation Indicator turns mostly red.**

9. If targeting multiple locations, repeat the above steps for each site.

5. FINALIZING THE PROCEDURE

1. Verify Balloon Deflation

Before removing the device, confirm that the Balloon is fully deflated.

- The Inflation Indicator should be mostly red.

NOTE	If intervention is needed during withdrawal (e.g., biopsy, polypectomy), the Balloon can be re-inflated for stabilization
-------------	---

2. Straighten the endoscope's tip and withdraw the endoscope from the device carefully.

NOTE	If fluids accumulate between the endoscope and Sleeve , remove the endoscope carefully to minimize spillage.
-------------	---

CAUTION	Do not use additional tools to loosen or remove the device—this may damage the endoscope.
----------------	--

3. Proper Disposal Procedure

After use, the DiLumen EZ¹ device is contaminated.

- Handle and dispose of the device according to hospital protocols and local, state, and federal regulations.

Dispose of all components, including:

- Sleeve Assembly
- Inflation Handle
- Inflation Line

CAUTION	Do not re-use. Single-use only.
----------------	---------------------------------

6. TROUBLESHOOTING

Issue: Balloon Cannot be Deflated

1. Ensure the Pinch Clamp is open.
2. Inspect the Inflation Line for blockages:
 - Ensure the Inflation Line is not bent or kinked.
 - Verify the Inflation Line Connector is securely attached.
 - Try deflating again:
 - Turn the Selection Switch to deflate (-).
 - Squeeze the Inflation Bulb slowly until the Inflation Indicator turns mostly red.
3. If the Balloon remains inflated:
 - Disconnect the Inflation Line Connector and check if the Balloon deflates.
 - Gently adjust the device position to assist deflation.
4. If deflation still fails:
 - Insert a small-diameter endoscope alongside the primary endoscope to visualize the Balloon.
 - Use an endoscopic tool to puncture the Balloon.
 - Confirm Balloon deflation through the small-diameter endoscope.
 - Remove the small-diameter endoscope and any tools.
 - Carefully withdraw the primary endoscope and the device.

7. STORAGE

Store DiLumen EZ¹ at room temperature in a dry environment.

CAUTION	Do not use past the expiration date printed on the label
----------------	--

8. SYMBOL LEGEND

	Medical Device		Do not use if package is damaged or open and consult instructions for use
	Authorized Representative in the European Community / European Union		Do not re-use
	Quantity Symbol		Temperature Limit
	Federal (US) law restricts this device to sale by or on the order of a physician		Humidity Limitation
	Consult Instructions for use or consult electronic instructions for use		Keep dry
	Use-by date		Medical Device Manufacturer
	Catalog Number		Caution
	Lot Number / Batch Code		Non-sterile

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