


Lumendi DiLumen Specifications

Intended Use / Indications for Use	The Lumendi DiLumen is an accessory to an endoscope. The DiLumen dual balloon accessory is intended for use with any endoscope that has an outer diameter of 12.5 - 14.3 mm and a minimum working length of 1680 mm. The device is indicated to ensure complete positioning of an endoscope in the large intestine, and assist with optical visualization, diagnosis, and endoscopic treatment.
Sterile/non-sterile	Non-sterile
Single use/ reusable	Single Use
Balloon Material	Low durometer Polyurethane
Sleeve/Overtube Material	Pellethane 2363-80A with 2.5% Glycolube translucent
Inflation tube/Pushrod Assembly Material	Isoplast 2510 Pebax 7233 SA01 Med Pebax 4033 SA01 Med Pellethane 90A
Compatible Endoscope Tip Diameter	12.5 - 14.3 mm
Fore Balloon Outer Diameter	60 mm
Aft Balloon Outer Diameter	60 mm
Balloon Pressure	45 ± 12 mmHg
Balloon Relief Pressure	55 mmHg
Sleeve Outer Diameter	19.5 mm
Sleeve Inner Diameter	17 mm
Working Length	~1680 mm
Inflation Source	Manual inflation bulb
Shelf Life	6 months

DiLumen
Endolumenal Interventional Platform



DiLumen

Endolumenal Interventional Platform



253 Post Road West, Westport, CT 06880
Tel: (203) 436-2669 Fax: (203) 557-0459
Lumendi.com