

DISCOVER DYNAMIC RETRACTION


The world's first endolumenal retraction system that provides precise grasping, multi-directional lifting and controlled tension during endoscopic treatment.

DiLumen ⊕ DR

Dynamic Retraction System


DEFINING THE NEXT ERA in Intestinal Therapy

DiLumen +DR Dynamic Retraction System



DiLumen +DR Dynamic Retraction System consists of a controllable sliding balloon assembly with two suture loops of 2cm and 5cm diameters. Depending on polyp position, either loop can be fixed to mucosal tissue via a common hemostasis clip and pushed or pulled for improved tissue exposure.

DiLumen +DR
Dynamic Retraction System



.....

DiLumen +DR provides:

- Effective retraction of tissue
- Continuous tissue tension
- Visualization of submucosal tissue plane
- Precise dissection



253 Post Road West, Westport, CT 06880
Tel: (203) 463-2669 Fax: (203) 557-0459
Lumendi.com